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Sunil Hwang; Todd Mullis; **Andrew Lee***
IMCS, Inc - Irmo, South Carolina

integrated micro-chromatography systems

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	Quantification 定量分析	Identification 鉴定分析
Mass spectrometer	TSQEndure, Thermo Fisher	Q-Exactive Thermo Fisher
Liquid chromatography	Vanquish UPLC	Ultimate 3000 nano-UHPLC
Mobile phase A	0.1% formic acid (F.A.) in water	0.1% F.A. in water
Mobile phase B	0.1% F.A. in acetonitrile	0.1% F.A. in acetonitrile
LC gradient	5% - 35% B for 15 minutes	2% - 30% B for 170 minutes
Trap column	N/A	Acclaim PepMap 100 (C_{18}, 5 µm, 100 Å, 300 µm X 5mm)
Analytical column	Syncronis C_{18}, 100 X 2.1 mm, 1.7 µm	Acclaim PepMap RSLC (C_{18}, 2 µm, 100 Å, 75 µm X 15 cm)
Column oven temperature	40 °C	Room temperature

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REFERENCES

1. *Modular stop and go extraction tips with stacked disks for parallel and multidimensional Peptide fractionation in proteomics. Ishihama Y, Rappsilber J, Mann M.J Proteome Res. 2006 Apr;5(4):988-94.*
2. *Protocol for micro-purification, enrichment, pre-fractionation and storage of peptides for proteomics using StageTips. Rappsilber J, Mann M, Ishihama Y. Nat Protoc. 2007;2(8):1896-906.*
3. *Sample clean-up strategies for ESI mass spectrometry applications in bottom-up proteomics: Trends from 2012 to 2016. Tubaon RM, HaddadPR, Quirino JP. Proteomics. 2017 Mar 8. doi: 10.1002/pmic.201700011. [Epub ahead of print] Review.*
4. *Recent advances in sample preparation techniques for effective bioanalytical methods. Kole PL, Venkatesh G, Kotecha J, Sheshala R. Biomed Chromatogr. 2011 Jan;25(1-2):199-217. doi: 10.1002/bmc.1560. Epub 2010 Dec 10. Review.*

Abbreviations:

ACN: Acetonitrile; AUC: Area under the curve; C_{18}A1: C_{18} 100 Å resin; C_{18}A3: C_{18} 300 Å resin; C_{18}N: C_{18} spin column from vendor N; C_{18}T: C_{18} spin column from vendor T; RP: wettable polystyrene cross linked with divinylbenzene; F.A.: Formic acid; GCB: Graphitized carbon black; TFA: Trifluoroacetic acid.

Contact: L. Andrew Lee - lee@imcstips.com

Artwork and layout created by Carmen Adamson.

© 2017 IMCS, Inc. All rights reserved. IMCStips® is a trademark of Integrated Micro-Chromatography Systems, Inc.